

Where Convenience Meets Compliance

EZ-DOCK SMARTBOX® Euthanasia Station



EZ-DOCK-OP

The EZ-Dock combines the best features of all our systems. A built-in flow circuit and times are already preset for the specific cage type to ensure AVMA compliancy. Plug in, connect the provided gas line and the system is ready for use! Once the cage is placed in dock properly, the cycles begin and shuts off when time has elapsed. The unit is designed to be as simple to use as possible.

Features:

- Compliant with AVMA Guidelines
- Cycle starts automatically when cage is docked
- Choose between adult and neonates
- LED indicator lights show when cycle is running and when complete
- Automatically shuts off if cage is removed for safety
- Uses house or bottled CO₂ supply at 15-100 psi
- Super E-Z

SMARTBOX
AUTO CO₂ SYSTEM



Scan to
watch video

