



**EZ-AF9000** - shown with E-Z Heat Controller and Heated Surgery Beds

### **EZ-AF9000 Auto Flow System**

Precise, consistent, preset flow rate from each outlet, no variables...

The E-Z Anesthesia Auto Flow System offers highly efficient, exact gas flow that is easy to use.

- Ideal for work with multiple animals and unique configurations
- Precise, low flow gas saves money and protects personnel
- Flow rates are set specifically for mice and rats
- Five outputs supply anesthesia to the induction chamber and up to four breathers
- Preset flow rates eliminate the need for flow meters and manual adjustments
- Protocols are easily followed without calculations
- Flow rate is not impacted by the type or the number of accessories connected
- Minimal opportunity for operator error
- Each output is individually controlled with a simple on/off switch
- Pushbutton oxygen flush purges induction chamber before removing animals
- Carry handle facilitates safe, easy transport of the unit

### **EZ-AF9000-A Customized Auto Flow System**

- Outputs can be individually customized to specific flow rates
- Can be customized to accommodate guinea pigs and rabbits

#### **Specifications:**

Anesthetic: Isoflurane 0 - 5%

Carrier Gas: Oxygen or Air

Flow Rate: Chamber - 1 LPM each

Breathers - 0.5 LPM each

Number of Outputs: 5

**EZ-AF9000 Base System Includes:**

EZ-190F Auto-Flow Anesthesia Machine

(1) Oxygen or Air Regulator

(2) EZ-103A Microflex Breathers

(1) EZ-178 Sure-Seal Mouse/Rat Induction Chamber


(1) Case of ReFresh Charcoal Filters (8)

Also includes mouse/rat nosecones, nosecone plugs, plastic surgery bed, and all connecting tubes

**EZ-AF9000  
Base System**



**EZ-AF9000 Auto Flow System**



- A. Oxygen Tank (not included)
- B. EZ-320 Oxygen Regulator
- C. EZ-15021 Oxygen Hose
- D. EZ-155 Tec 3 Isoflurane Vaporizer
- E. EZ-190F Auto-Flow Vaporizer Machine
- F. EZ-260N Connection Hose
- G. EZ-178 Mouse/Rat Sure-Seal Induction Chamber
- H. EZ-103A Microflex Mouse/Rat Non Rebreathing Device
- I. Plastic Surgery Bed
- J. EZ-251A Filter Connection Hose
- K. EZ-258 ReFresh Charcoal Filter
- L. HB-163 Electronic Heated Chamber Warmer (Optional)
- M. HB-103 Electronic Heating Bed with EZ-103H Breather (Optional)